



INSTRUMENT CARE

Sterion[®] | SafeClean[™] | Oxboro[™] | ColorTrack[™]

Instrument Care

Instrument Tip Protectors

Instrument Tip Protectors are designed to protect valuable instruments from breaking and chipping during sterilization and handling, help eliminate reprocessing costs caused by torn pouches or wrap, and protect staff from sharp sticks during handling.

Osteotome Tip Protectors

Save Osteotome instrument edges from chipping with these tip protectors. Available tinted or clear and vented or non-vented to comply with user preferences. All types allow visual instrument inspection and are gas and steam permeable. Latex-free and single-use.



Clear
50 Per Bag


Part No. Vented	Part No. Non-Vented	mm
094031bbg	094021bbg	6.4 x 19
094032bbg	094022bbg	9.5 x 25
094033bbg	094023bbg	12.7 x 25
094034bbg	094024bbg	19 x 25.4
094035bbg	094025bbg	25 x 25
094036bbg	094026bbg	32 x 32
094037bbg	094027bbg	38 x 38
094038bbg	094028bbg	44.5 x 32
094040bbg	094030bbg	Assorted

Blue
50 Per Bag

Part No. Vented	Part No. Non-Vented	mm
094011bbg	094001bbg	6.4 x 19
094012bbg	094002bbg	9.5 x 25
094013bbg	094003bbg	12.7 x 25
094014bbg	094004bbg	19 x 25.4
094015bbg	094005bbg	25 x 25
094016bbg	094006bbg	32 x 32
094017bbg	094007bbg	38 x 38
094018bbg	094008bbg	44.5 x 32
094020bbg	094010bbg	Assorted

Standard Tip Protector Sizing Guide

Use this guide to determine the **standard tip protector** on the following page that will best fit your instrument.



Standard Instrument Tip Protectors

Standard Instrument Tip Protectors are available solid, tinted or clear and vented or non-vented. They are gas or steam permeable and can be used in flash or conventional autoclave processes. All are available in 100 per package. Latex-free and single-use.



1.6mm x 19mm

Solid Vented	091011bbg
Solid Non-Vented	091001bbg
Tinted Vented	093011bbg
Tinted Non-Vented	093001bbg
Clear Vented	092011bbg
Clear Non-Vented	092001bbg



2mm x 19mm

Solid Vented	091012bbg
Solid Non-Vented	091002bbg
Tinted Vented	093012bbg
Tinted Non-Vented	093002bbg
Clear Vented	092012bbg
Clear Non-Vented	092002bbg



2.8mm x 19mm

Solid Vented	091013bbg
Solid Non-Vented	091003bbg
Tinted Vented	093013bbg
Tinted Non-Vented	093003bbg
Clear Vented	092013bbg
Clear Non-Vented	092003bbg



3.2mm x 25.4mm

Solid Vented	091014bbg
Solid Non-Vented	091004bbg
Tinted Vented	093014bbg
Tinted Non-Vented	093004bbg
Clear Vented	092014bbg
Clear Non-Vented	092004bbg



5mm x 25mm

Solid Vented	091015bbg
Solid Non-Vented	091005bbg
Tinted Vented	093015bbg
Tinted Non-Vented	093005bbg
Clear Vented	092015bbg
Clear Non-Vented	092005bbg



10mm x 19mm

Solid Vented	091016bbg
Solid Non-Vented	091006bbg
Tinted Vented	093016bbg
Tinted Non-Vented	093006bbg
Clear Vented	092016bbg
Clear Non-Vented	092006bbg



2mm x 9mm x 25mm

Solid Vented	091017bbg
Solid Non-Vented	091007bbg
Tinted Vented	093017bbg
Tinted Non-Vented	093007bbg
Clear Vented	092017bbg
Clear Non-Vented	092007bbg



2mm x 16mm x 25mm

Solid Vented	091018bbg
Solid Non-Vented	091008bbg
Tinted Vented	093018bbg
Tinted Non-Vented	093008bbg
Clear Vented	092018bbg
Clear Non-Vented	092008bbg



3mm x 25mm x 25mm

Solid Vented	091019bbg
Solid Non-Vented	091009bbg
Tinted Vented	093019bbg
Tinted Non-Vented	093009bbg
Clear Vented	092019bbg
Clear Non-Vented	092009bbg

Assorted — Solid

Solid Vented	091020bbg
Solid Non-Vented	091010bbg

Assorted — Tinted

Tinted Vented	093020bbg
Tinted Non-Vented	093010bbg
Tinted Vented Micro	093000bbg

Assorted — Clear

Clear Vented	092020bbg
Clear Non-Vented	092010bbg

Instrument Care

DuoGuards™

Ensure that your instruments stay in the proper open position during processing with DuoGuards™. Available vented and non-vented, DuoGuards™ are tinted for easy viewing of instrument tips. Latex-free and single-use.



Tinted
50 Per Bag

Part No.	Part No.	Color	mm
Vented	Non-Vented		
093031bbg	093021bbg	White	6.4 x 15
093032bbg	093022bbg	Blue	9.5 x 20
093033bbg	093023bbg	Green	12.7 x 20
093034bbg	093024bbg	Red	19 x 25.5
093035bbg	093025bbg	Yellow	15 x 25
093037bbg	093027bbg	Orange	25 x 25
093040bbg	093030bbg	Assorted	Assorted

InnerVent™ Instrument Tip Protectors

InnerVent™ Instrument Tip Protectors offer the performance of vented tip protectors with the safety of non-vented tip protectors. The raised inner nubs allow improved flow through of steam and gas while preventing the tip protectors from “vacuuming” onto the instrument during drying. Latex-free and single-use.



Tinted



Solid/Radiopaque



Clear

Solid, Clear, Tinted
100 Per Bag

Part No.	Color	mm
091117bbg	Orange, Solid	2 x 9 x 25
091118bbg	Brown, Solid	2 x 16 x 25
091119bbg	Black, Solid	3 x 25 x 25
092117bbg	Clear	2 x 9 x 25
092118bbg	Clear	2 x 16 x 25
092119bbg	Clear	3 x 25 x 25
093117bbg	Orange, Tinted	2 x 9 x 25
093118bbg	Brown, Tinted	2 x 16 x 25
093119bbg	Black, Tinted	3 x 25 x 25

Low Temperature Web Guards

Another Aspen exclusive, these unique guards with a flow through, open weave are designed to be used exclusively with low heat, plasma sterilization processes, such as the Sterrad™* Sterilization System. They are great for protecting endoscopes and their accessories. Latex-free and single-use.



*Sterrad is a Trademark of Johnson & Johnson Medical, Inc.

Part No.	Color	mm	Quantity
099001bbg	Orange Mesh	10 x 76	24
099002bbg	Green Mesh	15 x 76	24
099003bbg	Blue Mesh	20 x 76	24
099004bbg	Red Mesh	40 x 76	24
099005bbg	Black Mesh	60 x 76	16
099006bbg	White Mesh	70 x 76	16
099007bbg	Assorted Mesh	Assorted	20

Gelpi Retractor Tip Protectors

Protect expensive Gelpi retractors with Aspen Surgical's exclusive Gelpi Retractor Tip Protectors.

Available in 3 sizes to fit the most commonly used Gelpi retractors, these tip protectors are gas and steam sterilization compatible. Latex-free and single-use.



Part No.	Color	mm	Quantity Per Bag
096109bbg	Clear	38 Length	55
096129bbg	Clear	51 Length	45
096139bbg	Clear	76 Length	40

Weitlaner Instrument Tip Protectors

The Weitlaner Retractor Tip Protector is specifically designed for use with up to a 4" Weitlaner retractor and offers easy one-handed application. Protect your expensive Weitlaner instruments and reduce reprocessing costs associated with compromised wrap. Latex-free and single-use.



Vented, Solid or Clear
50 Per Bag

Part No.	Size	Color
096409bbg	Small	Magenta
096410bbg	Small	Clear
096419bbg	Medium	Green
096420bbg	Medium	Clear

Laparoscopic MIS Tip Protectors

Venting of these specialty tip protectors provides extra flow through to Trocar probes and Veress needles during high level disinfection or sterilization. Latex-free and single-use.



Vented and Clear

Part No.	Applications	Quantity Per Bag
097011bbg	19mm, Use for Veress Needle	50
097012bbg	20mm 5mm, Use for Trocar Probe	25
097013bbg	43mm 5mm, Use for Trocar Probe	20
097014bbg	20mm 10mm, Use for Trocar Probe	25
097015bbg	57mm 10mm, Use for Trocar Probe	20

Endoscope Guards

Safeguard the eyepieces of delicate endoscopic equipment during storage with these soft, conformable guards. Latex-free and single-use.



Blue, Non-Vented
6 Per Bag

Part No.	Description
095016bbg	Small Eye Piece

Round Tip Protectors

Designed to protect valuable instruments ranging in sizes from the smallest Veress Needle to the largest endoscope during handling and storage. Tinted for easy visual inspection and vented for improved steam, gas, or disinfectant flow. Latex-free and single-use.



Vented and Tinted

Part No.	Part No.	Color	Fits Tip Sizes mm	Quantity Per Bag
38.1mm	50.8mm			
096001bbg	096021bbg	White	2-3	70
096002bbg	096022bbg	Blue	3-4	70
096003bbg	096023bbg	Green	4-5	70
096004bbg	096024bbg	Red	5-6	70
096005bbg	096025bbg	Yellow	6-7	70
096006bbg	096026bbg	Maroon	7-8	70
096007bbg	096027bbg	Orange	8-9	70
096008bbg	096028bbg	Brown	9-10	70
096009bbg	096029bbg	Black	10-11	70
096013bbg	096033bbg	Green	15-16	50
096015bbg	096035bbg	Yellow	19-20	50
096017bbg		Orange	21-22	50

Instrument Care

Instrument Identification Tape

Tape makes organizing and tracking valuable instruments easy! Tape is flexible, adheres securely to instruments, and will not chip, flake, peel, or crack — and with 75+ colors, symbols, and character options, the color-coding options are nearly endless. If a change is needed, simply remove the ColorTrack I.D. Tape and replace with your new color choice. Latex-free.

Sheet Tape

Each 8½" x 11" sheet has 374 total inches of tape and contains:

- 44 each of 2 inch strips
- 44 each of 1½ inch strips
- 220 each of 1 inch strips

Solid Color

¼" Height / Various widths

Part No.	Description
152001eea	White
152002eea	Royal Blue
152003eea	Kelly Green
152004eea	Red
152005eea	Yellow
152006eea	Purple
152007eea	Orange
152008eea	Brown
152009eea	Black
152010eea	Gray
152011eea	Lavender
152012eea	Beige
152013eea	Light Blue
152014eea	Lime Green
152015eea	Pink
152017eea	Plum
152018eea	Mint Green

Dots

¼" Height / Various widths

Part No.	Description
158040eea	Dots, Fuchsia
158041eea	Dots, Medical Blue
158042eea	Dots, Emerald Green
158043eea	Dots, Violet

Stars

¼" Height / Various widths

Part No.	Description
158030eea	Stars, Fuchsia
158031eea	Stars, Medical Blue
158032eea	Stars, Emerald Green
158033eea	Stars, Violet

White Stripe

¼" Height / Various widths

Part No.	Description
154002eea	White/Royal Blue
154003eea	White/Kelly Green
154004eea	White/Red
154005eea	White/Yellow
154007eea	White/Orange
154008eea	White/Brown
154011eea	White/Lavender
154013eea	White/Light Blue
154014eea	White/Lime Green
154015eea	White/Pink
154017eea	White/Plum

Diagonal Stripes

¼" Height / Various widths

Part No.	Description
159001eea	Diagonal, Purple/Black
159002eea	Diagonal, Orange/Black
159003eea	Diagonal, Blue/Black
159004eea	Diagonal, Rubine/Black
159005eea	Diagonal, Yellow/Black
159006eea	Diagonal, Green/Black
159007eea	Diagonal, White/Black

Tri-Color Stripe

⅜" Height / Various widths

Part No.	Description
155051eea	Red/White/Blue
155052eea	Red/White/Black
155053eea	Red/Yellow/Black
155054eea	Red/Black/Green

Black Stripe

¼" Height / Various widths

Part No.	Description
153001eea	Black/White
153003eea	Black/Kelly Green
153004eea	Black/Red
153005eea	Black/Yellow
153007eea	Black/Orange
153010eea	Black/Gray
153011eea	Black/Lavender
153012eea	Black/Beige
153013eea	Black/Light Blue
153014eea	Black/Lime Green
153015eea	Black/Pink
153017eea	Black/Plum

Triangles

¼" Height / Various widths

Part No.	Description
158050eea	Triangles, Fuchsia
158051eea	Triangles, Medical Blue
158052eea	Triangles, Emerald Green
158053eea	Triangles, Violet

Numbers

¼" Height / Various widths

Each sheet contains numbers 0–8

Part No.	Description
11111	Numbers, Fuchsia
22222	Numbers, Medical Blue
33333	Numbers, Emerald Green
44444	Numbers, Violet

Assorted Symbols ¼" Height / Various widths

Part No.	Description
158060eea	Red Hearts on White
158070eea	Smiley Face, Yellow/Black
158080eea	Alphabet, Fuschia/Wt, A-E
158061eea	White Hearts on Red

4-Color ¼" Height / Various widths

Part No.	Description
152031eea	Red/Yellow/Purple/Orange
152032eea	Blue/White/Black/Green

Specialty Stripe ¼" Height / Various widths

Part No.	Description
156003eea	Blue/Green
156004eea	Blue/Red
157002eea	Yellow/Blue
157003eea	Yellow/Green
157004eea	Yellow/Red

Code to Specialty ¼" Height / Various widths

Part No.	Description
ENT 158090eea	Mint/White ENT
ER ER 158091eea	Purple/White ER
GEN 158092eea	Royal Blue/White GEN
OB/GYN 158093eea	Lt Blue/White OB/GYN
ORTHO 158094eea	Kelly Green/White ORTHO
NEURO 158095eea	Orange/White NEURO
CARDIAC 158096eea	Red/White CARDIAC

Roll Tape Roll Tape ¼" x 300"

Part No.	Description
151001eea	White
151002eea	Royal Blue
151003eea	Kelly Green
151004eea	Red
151005eea	Yellow
151006eea	Purple
151007eea	Orange
151008eea	Brown
151009eea	Black

Roll Tape ⅛" x 250"

Part No.	Description
151011eea	White
151012eea	Royal Blue
151013eea	Kelly Green
151014eea	Red
151015eea	Yellow
151016eea	Purple
151017eea	Orange
151018eea	Brown
151019eea	Black



Application Instructions:

For secure placement through multiple processing cycles, ColorTrack I.D. Tape should be wrapped 1½ times around the Instrument and excess length should be trimmed. Instrument sterility under ColorTrack I.D. Tape is assured after processing with properly applied tape.*

* Test Results Available On Request

Size and color of tapes shown may vary from actual products.

Same Day Tape Remover

An easy, instrument friendly solution to remove ColorTrack Tape. Eliminates the risk of cuts and sticks common with scalpel or utility knife tape removal. Neutral scent makes Same Day Tape Remover comfortable to use.

Part No.	Description
151066eea	1 Gallon Container

Easily Removes
Our Tape!



Instrument Care

Write-On ID Labels

Ideal for customization of trays, storage areas, cabinets or carts, these color-coded labels are formatted on 8½"x 11" sheets with three different label sizes on each sheet. For best results use with permanent markers only. Latex-free.



Part No.	Description
161002eea	Blue Border
161003eea	Green Border
161004eea	Red Border
161005eea	Yellow Border

Color-Coded Identification Locking Tags

These tamper evident locking tags can be customized to a variety of uses such as identifying instruments, carts, trays, containers or cabinets. Available pre-printed with "Repair", "Sharpen", or "Good" on one side and writable on the other or blank on one side and writable on the other. Sold in bags of 100. Latex-free.



Part No.	Description
251001bbg	Red, Imprint: Repair, Write-On
251002bbg	Yellow, Imprint: Sharpen, Write-On
251003bbg	Green, Imprint: Good, Write-On
251004bbg	Blue, Write-On
251005bbg	White, Write-On
251006bbg	Green, Write-On
251007bbg	Black, Write-On
251008bbg	Purple, Write-On
251009bbg	Orange, Write-On
251010bbg	Assorted Colors, Write-On

Utility Markers

Our utility markers provide a permanent method of identifying instruments, specimens, and other objects. The markers are waterproof, non-smearing, and can be used on most surfaces. Utility markers are available with or without nine black and white medication labels or twenty color medication labels and come in regular or taper tip depending on your need. Latex-free.




Part No.	Description	Quantity Per Box
2850	Marker and 9 Label Set — Sterile	50
2850BN	Marker and 9 Label Set — Non-Sterile	100
2851	Marker — Sterile	50
2851BN	Marker — Non-Sterile	100



Part No.	Description	Quantity Per Box
2830	Marker and 9 Label Set — Sterile	50
2830BN	Marker and 9 Label Set — Non-Sterile	100
2831	Marker — Sterile	50
2831BN	Marker — Non-Sterile	100
ML5832	Marker and 20 Medication Label Set — Sterile	50

Corner Protectors

Protect your instrument trays from rips, holes, and tears in sterilization wrap due to sharp tray corners, improper handling, or heavy surgical trays. Aspen Surgical's Corner Protectors are designed to create a flexible barrier between the tray's sharp corners and the sterilization wrap. They can be used with steam, gas, ozone, and low temperature plasma sterilization and are reusable. Latex-free.



Part No.	Description	Quantity Per Bag
CP-80NS	Corner Protectors	80

SafeClean™ Instrument and Channel Cleaning Brushes are designed to clean surgical/medical instruments and endoscopic devices. The complete line of SafeClean™ Cleaning Brushes are strong enough to effectively remove surgical residue, yet gentle enough to use on the most delicate instruments and devices. Latex-free.

General Purpose Brushes

General Purpose Brushes can be used throughout all of your instrument and device reprocessing areas to clean uniquely shaped and larger bore surgical, medical, and endoscopic instruments and devices. Latex-free.



Part No.	Brush Head Diameter	Brush Total Length	Bristle Length	Quantity Per Bag
241007bbg	1mm (0.039")	18"	0.5"	2
241010bbg	1.6mm (0.063")	18"	0.5"	2
241015bbg	1.8mm (0.071")	18"	0.5"	2
241019bbg	2mm (0.079")	12"	1"	2
241020bbg	2.2mm (0.086")	18"	0.5"	2
241024bbg	2.5mm (0.098")	24"	2.5"	2
241025bbg	2.7mm (0.106")	18"	1"	2
241029bbg	3mm (0.118")	12"	1"	2
241030bbg	3.2mm (0.126")	18"	2"	2
241032bbg	4mm (0.158")	12"	1"	2
241033bbg	5mm (0.197")	12"	2.5"	2
241034bbg	5mm (0.197")	24"	2.5"	2
241035bbg	5mm (0.197")	16"	2.5"	2
241040bbg	6mm (0.236")	16"	2.5"	2
241041bbg	6mm (0.236")	24"	2.5"	2
241042bbg	7mm (0.275")	16"	2.5"	2
241043bbg	7mm (0.275")	24"	2.5"	2
241044bbg	10mm (0.393")	16"	2.5"	2
241045bbg	11mm (0.433")	16"	2.5"	2
241046bbg	11mm (0.433")	24"	2.5"	2
241047bbg	10mm (0.393")	24"	2.5"	2
241050bbg	12mm (0.472")	16"	2.5"	2
241055bbg	15mm (0.590")	16"	2.5"	2
241060bbg	20mm (0.787")	16"	2.5"	2

Kits

Make ordering easy with all of the brushes you need in one kit!

Small Cleaning Brush Kit

This kit includes various small general purpose cleaning brushes. Latex-free.

Part No.

241100bbg — 1-Each

241010bbg

241015bbg

241020bbg

241025bbg

241030bbg

General Cleaning Brush Kit

This kit contains an extensive selection of general purpose brushes, as well as toothbrush and fan tip brushes. Latex-free.

Part No.

241600bbg — 3-Each

241025bbg

241030bbg

241035bbg

241040bbg

241041bbg

241045bbg

241046bbg

241050bbg

241055bbg

241060bbg

242001bbg

242002bbg

242003bbg

243001bbg

243004bbg

Large Cleaning Brush Kit

This kit includes various large general purpose cleaning brushes. Latex-free.

Part No.

241200bbg — 1-Each

241035bbg

241040bbg

241045bbg

241050bbg

241055bbg

241060bbg

"Best Sellers" Cleaning Brush Kit

This kit includes a selection of our most popular general purpose cleaning brushes. Latex-free.

Part No.

241300bbg — 1-Each

241010bbg

241015bbg

241020bbg

241025bbg

241030bbg

241035bbg

241040bbg

241045bbg

241050bbg

241055bbg

241060bbg

Micro/Ophthalmic Brush Kit

This kit contains a selection of our smallest general purpose cleaning brushes. Latex-free.

Part No.

241700bbg — 3-Each

241007bbg

241010bbg

241015bbg

241020bbg

241025bbg

241030bbg

Instrument Care

Suction Tube Cleaning Brushes

Suction Tube Cleaning Brushes are designed to fit into and effectively clean the smallest bore devices and instruments used in surgery, MIS, and therapy procedures today. Latex-free.



Part No.	Brush Head Diameter	Brush Total Length	Bristle Length	Quantity Per Bag
241003bbg	3 French (1mm, 0.039")	6"	0.5"	2
241005bbg	5 French (1.67mm, 0.066")	6"	0.5"	2
241006bbg	6 French (2mm, 0.079")	8"	0.5"	2
241008bbg	7 French (2.3mm, 0.092")	8"	0.5"	2
241009bbg	8 French (2.7mm, 0.105")	8"	0.5"	2
241011bbg	9 French (3mm, 0.118")	8"	0.5"	2
241012bbg	10 French (3.3mm, 0.131")	8"	0.5"	2
241013bbg	11 French (3.7mm, 0.144")	8"	0.5"	2
241014bbg	12 French (4mm, 0.158")	8"	0.5"	2

Fan-Tip Brush

Fan-Tip Brushes feature bristles that wrap around the end of the brush for superior cleaning of bottles, jars, and other large instruments. Latex-free.



Part No.	Brush Head Diameter	Brush Total Length	Bristle Length	Quantity Per Bag
242001bbg	4.8mm-12.7mm Tapered (0.189"-0.5")	7"	3"	2
242002bbg	13mm (0.512")	6"	3"	2
242003bbg	16mm (0.623")	5"	2.5"	2
242046bbg	9.6mm (0.384")	24"	2.5"	2

Tracheostomy Brushes

Use the right brush for the right job. Our Tracheostomy Brushes are uniquely designed to clean tubes. Latex-free.



Part No.	Brush Head Diameter	Brush Total Length	Bristle Length	Quantity Per Bag
247001bbg	Fits Trach Sizes 00-5FR 9mm (0.354")	4.75"	1.85"	2
247002bbg	Fits Trach Sizes 6-12FR 13mm (0.512")	8"	3.66"	2

Toothbrush Style

Aspen Surgical offers a variety of bristle types, designed to clean all types of instruments and stains including nylon, brass, and steel. Latex-free.



Part No.	Brush Head Diameter	Brush Total Length	Bristle Length	Quantity Per Bag
243001bbg	13mm (0.512"), Nylon	7.25"	1.5"	2
243002bbg	13mm (0.512"), Brass	7.25"	1.5"	2
243003bbg	13mm (0.512"), Steel	7.25"	1.5"	2

Double-Ended Nylon Cleaning Brushes

Our Double-Ended Cleaning Brush features one end with an extra large brush head for general cleaning and one end with a single row of bristles for more intricate cleaning. Latex-free.




Available
Sterile and
Non-Sterile

Part No.	Brush Head Diameter	Brush Total Length	Bristle Length	Quantity
243004bbg	Non-Sterile Short End, 2mm (0.079") Long End, 8mm (0.315")	7"	0.867" 1.38"	2/Bag
243014pbx	Sterile Short End, 2mm (0.079") Long End, 8mm (0.315")	7"	0.867" 1.38"	1/Pack, 12/Packs/Box

Reusable Channel Cleaning Brushes


Channel Brushes with tapered brush heads, elevator sheath, and comfortable round, plastic handles are available in a variety of lengths and materials to fit the range of endoscopic devices. Latex-free.



Part No.	Description	Total Length	Sheath Diameter	Quantity Per Bag
248107bbg	3–5mm	70cm	1.6mm	1
248115bbg	3–5mm	115cm	1.6mm	1
248160bbg	3–5mm	160cm	1.6mm	1
248180bbg	3–5mm	230cm	1.8mm	1
248182bbg	3–5mm — Double Ended	230cm	1.8mm	1
248185bbg	3–5mm — Teflon Coated	230cm	1.8mm	1
248190bbg	3–5mm	120cm	1.8mm	1

Disposable Scope Channel Cleaning Brushes

Aspen Surgical’s new Disposable Scope Channel Cleaning Brushes have brush heads that range from to 2.5mm – 5mm and will fit scopes as small as 1.1mm. Our disposable scope channel cleaning brushes feature more than double the bristle length than competing low-cost brushes. They have a strong, sturdy handle making it easy to pull through the scope, while still cleaning the scope effectively. Single-use and latex-free.




Disposable Scope Channel Brushes

Part No.	Description	Total Length	Sheath Diameter	Quantity Per Box
248102bbg	2.5mm	230cm	1mm	50
248103bbg	3mm	230cm	1.55mm	50
248104bbg	5mm — Double Ended	230cm	1.85mm	50
248105bbg	5mm	230cm	1.85mm	50

Hand/Nail Brushes

Our Hand/Nail brush is used to clean dirt on hands and nails. It can also be used for general instrument cleaning. Latex-free.



Part No.	Brush Head Diameter	Quantity Per Bag
244001bbg	35mm x 93mm (1.375" x 1.75")	2

Instrument Brush Storage

Convenient storage with our acrylic racks helps organize your cleaning brushes and help ensure you will always have the right brush for the job. Latex-free.

Part No.	Description	Quantity
246001bbg	Brush Rack, Medium (without drip tray)	1
246003bbg	Brush Rack, Large (with drip tray)	1



In addition to Instrument Care, Aspen Surgical carries a full line O.R. Accessory, Woundcare, Patient and Staff Safety, and Ultracell® Ophthalmic and Otolaryngology products. Please call us to request one of the catalogs above, or to share your new product ideas!



6945 Southbelt Dr SE
Caledonia, MI 49316 USA
Phone: 616.698.7100
888.364.7004
Fax: 616.698.9281
Web: aspensurgical.com