

OPERATING ROOM ACCESSORIES

Dr. Fog® | Sterion® | RICHARD-ALLAN | WRITESITE® | Securline®

O.R. Accessories

Surgical Skin Markers

Our surgical marking pens provide an accurate and reliable way to mark the incision site prior to surgery. Pens contain prep-resistant gentian violet ink and are available with or without ruler and blank labels. The only exception to this is our Taper Tip All-Skin Marker which contains a deeper blue ink for demarcation on darker skin types. A variety of barrel and tip options are available depending on your need and preference. Single-use and latex-free.

Classic Fine Tip Skin Marker



Part No.	Description	Qty
2530	Marker, Ruler, and Label Set — Sterile	12/Box
2531	Marker — Sterile	12/Box
2630	Marker, Ruler, and Label Set — Sterile	50/Box
2631	Marker — Sterile	50/Box

Classic Regular Tip Skin Marker



Part No.	Description	Qty
2550	Marker, Ruler, and Label Set — Sterile	12/Box
2551	Marker — Sterile	12/Box
2650	Marker, Ruler, and Label Set — Sterile	50/Box
2651	Marker — Sterile	50/Box

Modern Regular Tip Skin Marker



Part No.	Description	Qty
2712	Marker, Ruler, and Label Set — Sterile	12/Box
2750	Marker, Ruler, and Label Set — Sterile	50/Box
2751	Marker — Sterile	50/Box

Taper Tip Skin Marker



Part No.	Description	Qty
2730	Marker, Ruler, and Label Set — Sterile	50/Box
2731	Marker — Sterile	50/Box

Taper Tip All-Skin Marker



Part No.	Description	Qty
2756	Marker, Ruler, and Label Set — Sterile	50/Box
2756BN	Marker — Non-Sterile	500/Case

Securline® Taper Tip Skin Marker



Part No.	Description	Qty
1000-00-PDG	Marker and Ruler Set — Sterile	10/Box, 10 Boxes/Case
1002-00-PDG	Marker and Ruler Set — Sterile	100/Box
100L-00-PDG	Marker, Ruler, and Label Set — Sterile	100/Box
1003-00-PDG	Marker and Ruler Set — Non-Sterile	100/Box

WRITESITE® Preoperative Marking Products

Our non-sterile markers or stamp pad with gentian violet ink are designed to be used in your preoperative patient surgical site marking protocol. The newest addition to this line, the WRITESITE® Jr., generates less waste for this single-use item. Both markers are offered in our regular tip. Single-use and latex-free.

WRITESITE® Preoperative Surgical Site Marker



Part No.	Description	Qty
2753BN	Marker — Non-Sterile	100/Box

WRITESITE® Jr. Small Preoperative Skin Marker



Part No.	Description	Qty
2900BN	Marker — Non-Sterile	100/Box

WRITESITE® Preoperative Surgical Site “YES” Stamp



Part No.	Description	Qty
WS-001NS	“Yes” Stamp — Non-Sterile	75/Box

Utility Markers

Our utility markers provide a permanent method of identifying instruments, specimens, and other objects. The markers have black ink, are waterproof, non-smearing, and can be used on most surfaces. Utility markers are available with or without nine black and white or twenty color medication labels and come in regular or taper tip depending on your need. We also offer the Multi-Ink marker, which has gentian violet skin ink on one end and black utility ink on the other. Latex-free.

Regular Tip Utility Marker



Part No.	Description	Qty
2850	Marker and 9 Label Set — Sterile	50/Box
2850BN	Marker and 9 Label Set — Non-Sterile	100/Box
2851	Marker — Sterile	50/Box
2851BN	Marker — Non-Sterile	100/Box

Taper Tip Utility Marker



Part No.	Description	Qty
2830	Marker and 9 Label Set — Sterile	50/Box
2830BN	Marker and 9 Label Set — Non-Sterile	100/Box
2831	Marker — Sterile	50/Box
2831BN	Marker — Non-Sterile	100/Box
ML5832	Marker and 20 Medication Label Set — Sterile	50/Box

Multi-Ink Marker



Part No.	Description	Qty
2754	Marker and 9 Label Set — Sterile	12/Box
2754BN	Marker and 9 Label Set — Non-Sterile	100/Box
2755	Marker — Sterile	12/Box
2755BN	Marker — Non-Sterile	100/Box

Select markers are available with Blank Labels, Rulers, or Medication Labels as noted.



Blank Labels



Rulers



20 Color Medication Labels



9 Black and White Medication Labels

Securline® Lab Markers

Our lab markers are versatile and can be used to make permanent markings on plastic, metal, and porcelain labware. The unique ink has the ability to write on labware even when the surface is damp. These markers are available in both fine and extra-fine tip as well as in black, blue, green, or red ink. Latex-free.

Securline® Extra-Fine Tip Lab Marker



Part No.	Description	Qty
1300-13-PDC	Blue Marker — Non-Sterile	10/Box
1300-15-PDC	Green Marker — Non-Sterile	10/Box
1300-16-PDC	Red Marker — Non-Sterile	10/Box
1300-20-PDC	Black Marker — Non-Sterile	10/Box

Securline® Fine Tip Lab Marker



Part No.	Description	Qty
1400-13-PDC	Blue Marker — Non-Sterile	10/Box
1400-15-PDC	Green Marker — Non-Sterile	10/Box
1400-16-PDC	Red Marker — Non-Sterile	10/Box
1400-20-PDC	Black Marker — Non-Sterile	10/Box

Securline® Superfrost Marker

The Superfrost Marker has permanent ink that resists common solvents used in histology and cytology. This pen has an extra-fine tip and is available in black or red ink. Latex-free.

Securline® Extra-Fine Tip Superfrost Marker



Part No.	Description	Qty
1450-16-PDC	Red Marker — Non-Sterile	10/Box
1450-20-PDC	Black Marker — Non-Sterile	10/Box

O.R. Accessories

Corneal Marking Pad

Aspen Surgical's Corneal Marking Pad can be used with many refractive/corneal marking instruments. This pre-inked, ready-to-use gentian violet marking pad saves prep time and allows you to focus on your patient. Sterile, single-use, and latex-free.



Part No.	Description	Quantity Per Box
40490	Corneal Marking Pad	10

Endoscopic Kittner

Our blunt dissecting instrument is 38cm long with a 5mm stainless steel shaft. It is available in both single and dual tip. The radiopaque tip is non-absorbent. The brushed finish allows for safe use during laser surgery. Sterile, single-use, and latex-free.



Part No.	Description	Quantity Per Box
KT-9101	Single Tip	20
KT-9105	Dual Tip	20

Dr. Fog® Original Anti-Fog


DF-3100, the original anti-fog solution, comes with our easy squeeze bottle and sponge kit. We also offer the DF-3120 Ready-To-Use Sponge which is dry to touch. The sponge is alcohol-free. Both configurations are sterile and latex-free.



Part No.	Description	Quantity Per Box
DF-3100	Original Solution with Sponge	20
DF-3120	Alcohol-Free Dry Treated Sponge	40

Cautery Tip Cleaner

Cautery Tip Cleaners assist in the cleaning and removing of material from cautery tips during surgical procedures. Our tip cleaner consists of a 2" x 2" foam pad with abrasive material on one side and strong adhesive on the other for attachment to both disposable and non-disposable drapes. It will clean both mono-polar and bi-polar probes, is radiopaque, latex-free, and intended for single-use.



Part No.	Quantity Per Box
032001pbx	35
TC-100	100
TC-100b	240

Vessel Loops

Vessel loops are used in occluding, retracting, and identifying arteries, veins, tendons and nerves in surgical procedures. Our vessel loops are extruded from medical-grade silicone, with the exception of the Richard-Allan vessel loops, which are made of high performance, medical grade elastomer. All vessel loops are latex-free, radiopaque, and intended for single-use.



2 Loops / Pkg; 10 Pkgs / Box Double Pouched

Part No.	Color	mm Width x Length	Inches Width x Length
011001pbx	Red Mini	0.83 x 406	0.03 x 16
011002pbx	Blue Mini	0.83 x 406	0.03 x 16
011003pbx	Yellow Mini	0.83 x 406	0.03 x 16
011004pbx	White Mini	0.83 x 406	0.03 x 16
011011pbx	Red Maxi	1.0 x 406	0.04 x 16
011012pbx	Blue Maxi	1.0 x 406	0.04 x 16
011013pbx	Yellow Maxi	1.0 x 406	0.04 x 16
011014pbx	White Maxi	1.0 x 406	0.04 x 16
011021pbx	White Supermaxi	1.3 x 559	0.05 x 22
011022pbx	Blue Supermaxi	1.3 x 559	0.05 x 22
011032pbx	Blue Micro	0.41 x 559	0.019 x 22
011301pbx	Red Mini	0.83 x 762	0.03 x 30
011302pbx	Blue Mini	0.83 x 762	0.03 x 30
011303pbx	Yellow Mini	0.83 x 762	0.03 x 30
011304pbx	White Mini	0.83 x 762	0.03 x 30
011311pbx	Red Maxi	1.0 x 762	0.04 x 30
011312pbx	Blue Maxi	1.0 x 762	0.04 x 30
011313pbx	Yellow Maxi	1.0 x 762	0.04 x 30
011314pbx	White Maxi	1.0 x 762	0.04 x 30

Vessel Loops continued

2 Loops / Pkg; 10 Pkgs / Box Carded and Pouched

Part No.	Color	mm Width x Length	Inches Width x Length
016001pbx	Red Mini	0.83 x 457	0.03 x 18
016002pbx	Blue Mini	0.83 x 457	0.03 x 18
016003pbx	Yellow Mini	0.83 x 457	0.03 x 18
016004pbx	White Mini	0.83 x 457	0.03 x 18
016011pbx	Red Maxi	1.0 x 457	0.04 x 18
016012pbx	Blue Maxi	1.0 x 457	0.04 x 18
016013pbx	Yellow Maxi	1.0 x 457	0.04 x 18
016014pbx	White Maxi	1.0 x 457	0.04 x 18

RICHARD-ALLAN Vessel Loops



2 Loops / Pkg; 50 Pkgs / Box Banded and Pouched

Part No.	Color	mm Width x Length	Inches Width x Length
3901	Maxi Red	1.0 x 457	0.04 x 18
3902	Maxi Blue	1.0 x 457	0.04 x 18
3903	Maxi White	1.0 x 457	0.04 x 18
3904	Maxi Yellow	1.0 x 457	0.04 x 18
3905	Mini Red	.83 x 457	0.03 x 18
3906	Mini Blue	.83 x 457	0.03 x 18
3907	Mini White	.83 x 457	0.03 x 18
3908	Mini Yellow	.83 x 457	0.03 x 18

4 Loops / Pkg; 25 Pkgs / Box

Part No.	Color	mm Width x Length	Inches Width x Length
3910	Four Color Maxi Set	1.0 x 457	0.04 x 18
3912	Four Color Mini Set	0.83 x 457	0.03 x 18

Clamp Covers

Clamp Covers assist in the atraumatic clamping of delicate vessels during surgery. They are available in radiopaque silicone, non-radiopaque fabric, and radiopaque fabric. We offer a variety of sizes to provide outstanding convenience and flexibility of use. These are packaged sterile or non-sterile, intended for single-use, and latex-free.




Part No.	Color	Material	Radiopaque	Length mm	Length Inches	Quantity Per Box
Sterile						
061011pbx	Red	Silicone	Yes	2.5 x 130	0.10 x 5	10 Pkgs of 2
061012pbx	White	Silicone	Yes	5.1 x 130	0.20 x 5	10 Pkgs of 2
061013pbx	Blue	Silicone	Yes	7.6 x 130	0.30 x 5	10 Pkgs of 2
071001pbx	White	Fabric	No	3.2 x 130	0.125 x 5	10 Pkgs of 4
071002pbx	White	Fabric	No	9.5 x 130	0.375 x 5	10 Pkgs of 4
072001pbx	White	Fabric	Yes	3.2 x 150	0.125 x 6	10 Pkgs of 4
072002pbx	White	Fabric	Yes	9.5 x 150	0.375 x 6	10 Pkgs of 4
072004pbx	White	Fabric	Yes	6.3 x 150	0.25 x 6	10 Pkgs of 4

Non-Sterile						
062011bbg	Red	Silicone	Yes	2.5 x 130	0.10 x 5	20
062012bbg	White	Silicone	Yes	5.1 x 130	0.20 x 5	20
062013bbg	Blue	Silicone	Yes	7.6 x 130	0.30 x 5	20

Suture Aid Booties

Suture Aid Booties protect, tag, and help locate sutures. Our vinyl suture booties are available in two-, three-, or five-pair configurations and are positioned in an adhesive-backed foam block for quick and easy usage. They are radiopaque, latex-free, and are packaged sterile for single-use. Choose from a variety of sizes and colors as listed below.



Part No.	Description	Quantity Per Tray	Quantity Per Box
051003pbx	Standard Yellow	5 Pair	5 Trays
051011pbx	Standard Red	5 Pair	5 Trays
051012pbx	Standard Blue	5 Pair	5 Trays
051014pbx	Standard White	5 Pair	5 Trays
051015pbx	Standard Asst.	5 Pair	5 Trays
051004pbx	Large Pink	5 Pair	5 Trays
053003pbx	Standard Yellow	3 Pair	5 Trays
052004pbx	Large Pink	2 Pair	5 Trays
052003pbx	Standard Yellow	2 Pair	5 Trays

O.R. Accessories

Bulldog Vascular Clamps

Aspen surgical provides three varieties of Bulldog Vascular Clamps for use in occluding vessels, depending on desired clamping force, jaw opening, and jaw length. All are sterile, single-use, and latex-free.

Atraumatic Bulldog Clamp

The jaws of the Atraumatic Bulldog Clamp are coated with a soft material that gently grips the vessel reducing trauma to the patient, while still successfully occluding the vessel. The clamps have virtually the same clamping force of competitors while requiring less force to apply.

- Clamping Force – 1.048 lbs
- Jaw Opening – 6mm
- Jaw Length – 19mm



Part No.	Description	Quantity Per Box
031001pbx	Atraumatic Bulldog Clamp	5 Pair

Secure-Hold Bulldog Clamp

The Secure-Hold Bulldog Clamp utilizes a soft material on the jaws of the clamp to occlude the vessel during use. This clamp is designed much like the Atraumatic Clamp to reduce a patient's trauma, while still occluding the vessel.

- Clamping Force – 1.149 lbs
- Jaw Opening – 6mm
- Jaw Length – 19mm



Part No.	Description	Quantity Per Box
031003pbx	Secure-Hold Bulldog Clamp	5 Pair

Extra-Traction Bulldog Clamp

The Extra-Traction Bulldog Clamp has a serrated design to occlude larger and tougher vessels. The clamp has a larger jaw opening than the competitor, while still maintaining similar opening and clamping forces.

- Clamping Force – 1.556 lbs
- Jaw Opening – 10mm
- Jaw Length – 19mm



Part No.	Description	Quantity Per Box
031004PBX	Extra-Traction Bulldog Clamp	5 Pair

Syringe and Scalpel Handling System

The Syringe and Scalpel Handling System is used to provide safe and easy handling and storage for syringes and scalpels. This item consists of an adhesive-backed block that firmly grips and holds the syringe cap, syringe, or scalpel. It is available with slots for one, two, or three syringes or scalpels, is intended for single use, and is latex-free.



Part No.	Block Size	Width mm	Width Inches	Quantity Per Box
Sterile				
111001pbx	Single	31.8	1.25	5
111002pbx	Double	44.5	1.75	5
111003pbx	Triple	69.9	2.75	5
Non-Sterile				
111001bbg	Single	31.8	1.25	12
111002bbg	Double	44.5	1.75	12
111003bbg	Triple	69.9	2.75	12

Pin Covers

Our Pin Covers help to protect patients, staff, and clothing from the sharp ends of Steinmann Pins and K-Wires. These are made of thermoplastic and are available in three color-coded sizes to fit all eleven sizes of Steinmann Pins and K-Wires. They are latex-free, sterile, and packaged for single-use.

- Yellow (Small) – Packaged 1/Pouch to match applications that involve one exposed end.
- Blue (Medium) and Red (Large) – Packaged 2/Pouch to match applications that involve two exposed ends.



Part No.	Color	Fits This Size Range/Inches	Quantity Per Box
101001pbx	Yellow	0.028 – 0.062	1/ Pouch, 60 Pouches/Box
101002pbx	Blue	0.078 – 0.125	2/Pouch, 30 Pouches/Box
101003pbx	Red	0.141 – 0.188	2/Pouch, 30 Pouches/Box
101000pbx	Assorted	Assorted	50 Assorted Pouches/Box

Disposable Light Handles

We offer two systems for convenience: Rigid Light Handle or Flexible Cover. Eliminate costly reprocessing with these easy to use disposable Rigid Handles or Flexible Covers. All are sterile, latex-free, and intended for single-use.

Adaptors are available for most of the popular lights. Please call for further adapter information.

Rigid Handle




Figure A:
Illustrates our Rigid Light Handle. Our threaded design fits most adaptors made.

Flexible Cover




Figure B:
The Flexible Cover was designed using our unique push system allowing you to easily apply the cover. The pull tab makes for fast removal.

Part No.	Description	Quantity Per Box
LT-R01	Rigid Handle	24
LT-F01b	Flexible Cover (1 Each)	80
LT-F02	Flexible Cover (2 Each)	80
LT-F03	Flexible Cover (3 Each)	80

Snares

Our tonsil and nasal snares are already formed in pre-cut, pre-bent sizes, ready for immediate use. Sterile, single-use, and latex-free.



Part No.	Description	Quantity Per Box
2405	Nasal Snare Size 5	2/Pack, 50 Packs/Box
2406	Nasal Snare Size 6	2/Pack, 50 Packs/Box
2407	Tonsil Snare Size 7	2/Pack, 50 Packs/Box
2408	Tonsil Snare Size 8	2/Pack, 50 Packs/Box
2409	Tonsil Snare Size 9	2/Pack, 50 Packs/Box

Safety Pins – Sterile


Our safety pins are packaged sterile, latex-free, and intended for single-use.



Part No.	Description Inches	Quantity Per Box
2201	1	2/Pack, 50 Packs/Box
2202	1.5	2/Pack, 50 Packs/Box
2203	2	2/Pack, 50 Packs/Box

Other Aspen Surgical Product Offerings

In addition to O.R. Accessories, Aspen Surgical carries a full line of Patient and Staff Safety, Woundcare, Instrument Care, and Ultracell® Ophthalmic and Otolaryngology products. For further information, please request any of the catalogs below.



Aspen Surgical is committed to on-going product development. Please call us to share your new product ideas!



6945 Southbelt Dr SE
Caledonia, MI 49316 USA
Phone: 616.698.7100
888.364.7004
Fax: 616.698.9281
Web: aspensurgical.com