



UNIMAX
UNIMAX MEDICAL SYSTEMS INC.



UNIMAX SURGICAL TROCARS

www.unimaxmeds.com

UNIMAX SURGICAL TROCAR



Thoracic trocar



Hasson trocar



Auto locking trocar



Secondary trocar




Bladeless trocar



Visible trocar

UNIMAX VISIBLE TROCAR

- Provides safe entry into the peritoneal cavity
- Reduces trocar-related injuries
- Reduces scarring
- No fascial closure is required
- Less pain and herniation rate



Clear cannula for better visualization
Ribbed cannula for optimal attachment



Accommodated with universal seal for the entry of 5-12mm instruments



Clear tip allows direct visualization of the abdominal wall layers



Detachable clear cannula head enables easier specimen retrieval

Ergonomic handle design ensures ease of insertion

SURGICAL TROCAR

PRODUCT CODE	DESCRIPTION	PACKAGING
FSTD05070	Auto locking trocar 5mm/7cm	5 pcs/box
FSTD05100	Auto locking trocar 5mm/10cm	
FSTD10100	Auto locking trocar 10mm/10cm-detachable cannula cap	
FSTD12100	Auto locking trocar 12mm/10cm-detachable cannula cap	
FST610150	Auto locking trocar 10mm/15cm	5 pcs/box
FST612150	Auto locking trocar 12mm/15cm	
FHTD10100	Hasson trocar 10mm/10cm-detachable cannula cap	5 pcs/box
FHTD12100	Hasson trocar 12mm/10cm-detachable cannula cap	
FPTD05070	Bladeless trocar 5mm/7cm	5 pcs/box
FPTD05100	Bladeless trocar 5mm/10cm	
FPTD10100	Bladeless trocar 10mm/10cm-detachable cannula cap	
FPTD12100	Bladeless trocar 12mm/10cm-detachable cannula cap	
FPT805071	Secondary trocar 5mm/7cm- bladeless	5 pcs/box
FVTD10100	Visible trocar 10mm/10cm-detachable cannula cap	5 pcs/ box
FVTD12100	Visible trocar 12mm/10cm-detachable cannula cap	
FTT106070	Thoracic trocar 6mm/7cm	10 pcs/box
FTT111070	Thoracic trocar 11mm/7cm	
FTT113070	Thoracic trocar 13mm/7cm	
FTT115070	Thoracic trocar 15mm/7cm	



UNIMAX
UNIMAX MEDICAL SYSTEMS INC.



UNIMAX SUCTION IRRIGATION SYSTEM

www.unimaxmeds.com

UNIMAX SUCTION IRRIGATION SYSTEM



SUCTION IRRIGATION SYSTEM

- Color-coded buttons and tubings for easy identification
- Various probes meet all kinds of demands
(available in 28cm, 33cm, 45cm working length and 5mm, 10mm diameter)



ELECTRODES WITH SUCTION CHANNEL



L hook



Spatula



Button electrode



Needle electrode

SUCTION IRRIGATION SYSTEM

PRODUCT CODE	DESCRIPTION	PACKAGING
FSI128050	Suction irrigation set 5mm/28cm	10 sets/box
FSI133050	Suction irrigation set 5mm/33cm	
FSI145050	Suction irrigation set 5mm/45cm	
FSI133100	Suction irrigation set 10mm/33cm	
FSI533050	Suction irrigation set 5mm/33cm with button electrode	10 sets/box
FSI633050	Suction irrigation set 5mm/33cm with L hook	
FSI833050	Suction irrigation set 5mm/33cm with spatula	
FPB105280	Disposable probe for S/I set 5mm/28cm	10 pcs/box
FPB105330	Disposable probe for S/I set 5mm/33cm	
FPB105450	Disposable probe for S/I set 5mm/45cm	
FPB110330	Disposable probe for S/I set 10mm/33cm	
1PB405330	Reusable probe for S/I set 5mm/33cm	10 pcs/box
1PB405450	Reusable probe for S/I set 5mm/45cm	
1PB410330	Reusable probe for S/I set 10mm/33cm	
FSA205330	Monopolar L hook probe for S/I set 5mm/33cm	10 pcs/box
FSA305330	Monopolar spatula probe for S/I set 5mm/33cm	
FSA505330	Monopolar button electrode probe for S/I set 5mm/33cm	
FSAN05330	Monopolar needle probe for S/I set 5mm/33cm	



UNIMAX
UNIMAX MEDICAL SYSTEMS INC.



UNIMAX MONOPOLAR SCISSORS

www.unimaxmeds.com

UNIMAX MONOPOLAR SCISSORS



MONOPOLAR SCISSORS

- Stainless steel blades sustain sharpness throughout the procedure
- Precise & sharp blades for smooth cuttings

360° rotation wheel allows the shaft to rotate for 360 degrees



Extended PTFE insulation reduces the chance of thermal damage

- Ergonomic handle increases comfort during operation
- Extra smooth manipulation

MONOPOLAR SCISSORS

PRODUCT CODE	DESCRIPTION	PACKAGING
FMS533050	Monopolar scissors 5mm/33cm	5 pcs/box



UNIMAX MEDICAL SYSTEMS INC.

8F-2, NO.127 Lane 235, Pao Chiao Rd, Hsin Tien, Taipei, Taiwan
T/+886-2-8919-1698 F/+886-2-8919-1528 www.unimaxmeds.com



UNIMAX
UNIMAX MEDICAL SYSTEMS INC.



UNIMAX
LAPAROSCOPIC
GRASPING
INSTRUMENTS

www.unimaxmeds.com

UNIMAX LAPAROSCOPIC GRASPING INSTRUMENTS



LAPAROSCOPIC GRASPING INSTRUMENTS

Monopolar cable connector in the rear of the handle

360° rotation wheel allows 360° shaft rotation

Quick release ratchet design

Ergonomic handle increases comfort during operation



Duckbill

Grasper

Clinch

Rat tooth

Dissector

Babcock

Johan

GRASPING FORCEPS

PRODUCT CODE	PRODUCT DESCRIPTION	PACKAGING
FEG233050	Straight grasping forceps 5mm/33cm	5 pcs/box
FED233050	Maryland dissecting forceps 5mm/33cm	5 pcs/box
FEC233050	Clinch grasping forceps 5mm/33cm	5 pcs/box
FEB233050	Babcock grasping forceps 5mm/33cm	5 pcs/box
FRT233050	Rat tooth grasping forceps 5mm/33cm	5 pcs/box
FDG233050	Duckbill grasping forceps 5mm/33cm	5 pcs/box
FEJ233050	Johan grasping forceps 5mm/33cm (38mm long clip)	5 pcs/box
FEJ23305A	Johan grasping forceps 5mm/33cm (22mm short clip)	5 pcs/box



UNIMAX MEDICAL SYSTEMS INC.

8F-2, NO.127 Lane 235, Pao Chiao Rd, Hsin Tien, Taipei, Taiwan
T/+886-2-8919-1698 F/+886-2-8919-1528 www.unimaxmeds.com



UNIMAX
UNIMAX MEDICAL SYSTEMS INC.



UNIMAX
ENDO
POUCH SERIES

www.unimaxmeds.com

UNIMAX ENDO POUCH SERIES



ENDO POUCH SERIES



Endo pouch with introducer



Endo pouch with suture tie



Endo pouch with memory wire

BENEFITS OF DETACHABLE ENDO POUCH

Unimax detachable endo pouch is a patented device which provides a quicker and easier application of endo pouch. This unique design allows the pouch to be detached from the introducer and left within Patient body during procedure. Consequently, it increases the convenience of procedure and shortens the procedure time.



- Unique patented design ensures ease of use
- 10mm introducer provides ease of insertion
- Self-opening pouch design provides easy and safe specimen capture
- High strength PU pouch material assures no leakage
- Enhanced metal ring ensures large pouch opening and keeps the pouch in place when loading specimen

ENDO POUCH SERIES

PRODUCT CODE	DESCRIPTION	PACKAGING
FEP136080	Endo pouch with introducer 3"*6"	20 pcs/box
FEP157080	Endo pouch with introducer 5"*7"	20 pcs/box
FEP736080	Endo pouch with suture tie 3"*6"	20 pcs/box
FEP757080	Endo pouch with suture tie 5"*7"	20 pcs/box
FEP836080	Endo pouch with memory wire 3"*6"	10 pcs/box
FEP857080	Endo pouch with memory wire 5"*7"	10 pcs/box
FEP936100	Detachable endo pouch 3"*6"	5 pcs/box
FEP957100	Detachable endo pouch 5"*7"	5 pcs/box



UNIMAX MEDICAL SYSTEMS INC.

8F-2, NO.127 Lane 235, Pao Chiao Rd, Hsin Tien, Taipei, Taiwan
 T/+886-2-8919-1698 F/+886-2-8919-1528 www.unimaxmeds.com



UNIMAX
UNIMAX MEDICAL SYSTEMS INC.



UNIMAX BIPOLAR INSTRUMENTS

www.unimaxmeds.com

UNIMAX BIPOLAR INSTRUMENTS



BIPOLAR INSTRUMENT SERIES



Bipolar dissecting forceps



Bipolar cutting loop

BIPOLAR CUTTING FORCEPS



Special coated jaws and blade remarkably reduce tissue adhesion

Jaw-locking design ensures secure grasp of tissues



- Multi-functional instrument incorporates grasping, dissecting, transecting, and coagulating functions in one
- Fewer instrument exchanges promote efficiency and reduce procedure time
- Offers various cable options for different bipolar generators

BIPOLAR INSTRUMENTS

PRODUCT CODE	DESCRIPTION	PACKAGING
FBS133050	Bipolar cutting forceps 5mm/33cm (US adaptor)	2 pcs/box
FBS133100	Bipolar cutting forceps 10mm/33cm (US adaptor)	
FBS233050	Bipolar cutting forceps 5mm /33cm (EU adaptor, two holes)	
FBS233100	Bipolar cutting forceps 10mm/33cm (EU adaptor, two holes)	
FBS333050	Bipolar cutting forceps 5mm/33cm (valleylab adaptor, two pins)	
FBS433050	Bipolar cutting forceps 5mm/33cm (three-pin adaptor)	
FBD133050	Bipolar dissecting forceps 5mm/33cm (US adaptor)	5 pcs/box
FBD233050	Bipolar dissecting forceps 5mm/33cm (EU adaptor)	
FBL100000	Bipolar cutting loop 4mm/28cm (US adaptor)	10 pcs/box
FBL200000	Bipolar cutting loop 4mm/28cm (EU adaptor)	
FBL300000	Bipolar cutting loop 4mm/28cm (PK adaptor)	



UNIMAX MEDICAL SYSTEMS INC.

8F-2, NO.127 Lane 235, Pao Chiao Rd, Hsin Tien, Taipei, Taiwan
T/+886-2-8919-1698 F/+886-2-8919-1528 www.unimaxmeds.com



UNIMAX
UNIMAX MEDICAL SYSTEMS INC.



UNIMAX
INSUFFLATION
TUBING
SET

www.unimaxmeds.com

UNIMAX INSUFFLATION TUBING SET



INSUFFLATION TUBING SET



15mm connector
(Wolf / Storz)



22mm connector
(Linvatec)



Click-in adaptor
(Stryker)



O ring adaptor
(Smith Nephew)



Luer Lock



15/22mm connector

- Compatible with most equipments
- Filter made of membrane:0.02 micron high quality filter(99.9999% filtration efficiency)
- Adequate hardness of the tubing prevents kinking
- Medical grade material used
- Six different connectors for options

INSUFFLATION TUBING SET

PRODUCT CODE	DESCRIPTION	PACKAGING
FIF231500	Insufflation tubing set with luer lock	10 pcs/box
FIF231510	Insufflation tubing set with O ring adaptor	10 pcs/box
FIF231520	Insufflation tubing set with click-in adaptor	10 pcs/box
FIF231530	Insufflation tubing set with 15mm round connector	10 pcs/box
FIF231540	Insufflation tubing set with 22mm round connector	10 pcs/box
FIF231550	insufflation tubing set with 15/22 mm round connector	10 pcs/box



UNIMAX MEDICAL SYSTEMS INC.

8F-2, NO.127 Lane 235, Pao Chiao Rd, Hsin Tien, Taipei, Taiwan
T/+886-2-8919-1698 F/+886-2-8919-1528 www.unimaxmeds.com



UNIMAX
UNIMAX MEDICAL SYSTEMS INC.



UNIMAX VERESS **NEEDLE**

www.unimaxmeds.com

UNIMAX VERESS NEEDLE



VERESS NEEDLE



Sharp tip for easy insertion



Red visual needle position indicator

- Two sizes of needle available
- High flow stopcock
- Spring-loaded stylet provides audible "click" when passing through tissue layers

VERESS NEEDLE

PRODUCT CODE	DESCRIPTION	PACKAGING
FVN112000	Veress needle 120mm	25 pcs/box
FVN115000	Veress needle 150mm	25 pcs/box



UNIMAX MEDICAL SYSTEMS INC.

8F-2, NO.127 Lane 235, Pao Chiao Rd, Hsin Tien, Taipei, Taiwan
T/+886-2-8919-1698 F/+886-2-8919-1528 www.unimaxmeds.com



UNIMAX
UNIMAX MEDICAL SYSTEMS INC.



UNIMAX LAPAROSCOPIC **ELECTRODES**

www.unimaxmeds.com

UNIMAX LAPAROSCOPIC ELECTRODES



LAPAROSCOPIC ELECTRODES



- Enables a lasting and precise coagulation
- Three different tip designs: L hook, J hook, spatula

LAPAROSCOPIC ELECTRODES

PRODUCT CODE	DESCRIPTION	PACKAGING
FPL533050	Monopolar L hook probe 5mm/33cm	10 pcs/box
FPL533052	Monopolar L hook probe 5mm/33cm with cable	10 pcs/box
FPJ533050	Monopolar J hook probe 5mm/33cm	10 pcs/box
FPJ533052	Monopolar J hook probe 5mm/33cm with cable	10 pcs/box
FPS533050	Monopolar spatula probe 5mm/33cm	10 pcs/box
FPS533052	Monopolar spatula probe 5mm/33cm with cable	10 pcs/box




UNIMAX MEDICAL SYSTEMS INC.

8F-2, NO.127 Lane 235, Pao Chiao Rd, Hsin Tien, Taipei, Taiwan
 T/+886-2-8919-1698 F/+886-2-8919-1528 www.unimaxmeds.com



UNIMAX
UNIMAX MEDICAL SYSTEMS INC.




UNIMAX CAMERA SLEEVES

www.unimaxmeds.com

UNIMAX CAMERA SLEEVES



CAMERA SLEEVES



Closed camera sleeve



Standard camera sleeve



Folded camera sleeve

Unimax disposable camera sleeve comes in 3 different models to satisfy various needs.

- Offers effective barrier
- Allows quick and easy application
- Prolongs the life of camera head and cable

CAMERA SLEEVES

PRODUCT CODE	DESCRIPTION	PACKAGING
FCS125000	Standard camera sleeve	20 pcs/box
FCS225000	Closed camera sleeve	20 pcs/box
FCS525000	Folded camera sleeve	50 pcs/box



UNIMAX

UNIMAX MEDICAL SYSTEMS INC.



UNIMAX SUCTION **LINER** **SYSTEM**

www.unimaxmeds.com

UNIMAX SUCTION LINER SYSTEM



SUCTION LINER SYSTEM



Shut-off valve

- Protection devices: with bacteria filter and shut off valve for better protection of patients and hospital device
- Safe: completely sealed and compact endurance
- Versatile: it can be used individually in the wall plate or used in tandem in the mobile suction cart(2~4 units)
- User friendly: color coded lid to determine the connection of suction liners when used in mobile suction cart
- Environmentally friendly: polyethylene used

SUCTION LINER SYSTEM

PRODUCT CODE	DESCRIPTION	PACKAGING
SL-1500	Disposable suction liner 1500c.c. (white lid)	50 pcs/carton
SL-1500S	Disposable suction liner 1500c.c. with shutoff valve(blue lid)	50 pcs/carton
SL-2000	Disposable suction liner 2000c.c. (white lid)	50 pcs/carton
SL-2000S	Disposable suction liner 2000c.c. with shutoff valve(blue lid)	50 pcs/carton
SL-1500CH	Reusable outer canister 1500c.c.	40 pcs/carton
SL-2000CH	Reusable outer canister 2000c.c.	16 pcs/carton




UNIMAX MEDICAL SYSTEMS INC.

8F-2, NO.127 Lane 235, Pao Chiao Rd, Hsin Tien, Taipei, Taiwan
T/+886-2-8919-1698 F/+886-2-8919-1528 www.unimaxmeds.com



UNIMAX
UNIMAX MEDICAL SYSTEMS INC.



UNIMAX
YANKAUER
SUCTION
PROBE

www.unimaxmeds.com

UNIMAX YANKAUER SUCTION PROBE



YANKAUER SUCTION PROBE



- Available in four types: rose tip with vent, rose tip without vent, flat tip with vent and flat tip without vent
- Made by medical grade resin which meets up to the requirement of Class VI

YANKAUER SUCTION PROBE

PRODUCT CODE	DESCRIPTION	PACKAGING
FYS130040	Yankauer suction probe rose tip without vent.	50 pcs/box
FYS130050	Yankauer suction probe rose tip with vent.	50 pcs/box
FYS230140	Yankauer suction probe flat tip without vent.	50 pcs/box
FYS230150	Yankauer suction probe flat tip with vent.	50 pcs/box



UNIMAX MEDICAL SYSTEMS INC.

8F-2, NO.127 Lane 235, Pao Chiao Rd, Hsin Tien, Taipei, Taiwan
T/+886-2-8919-1698 F/+886-2-8919-1528 www.unimaxmeds.com



UNIMAX
UNIMAX MEDICAL SYSTEMS INC.




UNIMAX LLETZ ELECTRODES

www.unimaxmeds.com

UNIMAX LLETZ ELECTRODES



LLETZ ELECTRODES



- Comes in various types and different sizes
- Reusable performance with much friendly price
- Unique design of the ball electrode prevents the ball from falling off during procedure
- Special coating ensures no adhesion of the ball electrode

LLETZ ELECTRODES

PRODUCT CODE	DESCRIPTION	PACKAGING
FLE110050	Disposable loop electrode 10mm*5 mm	50 pcs/box
FLE110100	Disposable loop electrode 10mm*10 mm	
FLE115080	Disposable loop electrode 15mm*8 mm	
FLE115100	Disposable loop electrode 15mm*10 mm	
FLE115120	Disposable loop electrode 15mm*12 mm	
FLE115150	Disposable loop electrode 15mm*15 mm	
FLE120080	Disposable loop electrode 20mm*8 mm	
FLE120100	Disposable loop electrode 20mm*10 mm	
FLE120120	Disposable loop electrode 20mm*12 mm	
FLE120150	Disposable loop electrode 20mm*15 mm	
FLE120200	Disposable loop electrode 20mm*20 mm	
FLE125150	Disposable loop electrode 25mm*15 mm	
FLE210100	Disposable square electrode 10mm*10 mm	50 pcs/box
FEF100700	Disposable blade electrode 7cm	50 pcs/box
FEF101250	Disposable blade electrode 12.5cm	
FEF101500	Disposable blade electrode 15cm	
FEN100700	Disposable needle electrode 7cm	50 pcs/box
FEN101500	Disposable needle electrode 15cm	
FEN200700	Disposable needle electrode -45 degree angle 7cm	
FCB109120	Cone biopsy excisor 9mm*12mm	50 pcs/box
FCB112150	Cone biopsy excisor 12mm*15mm	
FCB112200	Cone biopsy excisor 12mm*20mm	
FCB115180	Cone biopsy excisor 15mm*18mm	



UNIMAX MEDICAL SYSTEMS INC.

8F-2, NO.127 Lane 235, Pao Chiao Rd, Hsin Tien, Taipei, Taiwan
 T/+886-2-8919-1698 F/+886-2-8919-1528 www.unimaxmeds.com



UNIMAX
UNIMAX MEDICAL SYSTEMS INC.



UNIMAX IRRIGATION PUMP

www.unimaxmeds.com

UNIMAX IRRIGATION PUMP



IRRIGATION PUMP

Locking pressure adjustment knob



on/off switch

- Available in single 1 liter and dual 1 liter pump system
- Noise free
- Easy to operate
- Reduces O.R. time

IRRIGATION PUMP

PRODUCT CODE	DESCRIPTION	PACKAGING
1PP110000	Irrigation pump (reusable)	1 unit/box
1PP210000	Twin Irrigation Pump with mounting plate and clamp (reusable)	1 unit/box



UNIMAX
UNIMAX MEDICAL SYSTEMS INC.



UNIMAX ANTI-FOG

www.unimaxmeds.com



ANTI-FOG

With long lasting fog free formula, Unimax anti fog solution provides clear visualization for laparoscopic procedures.



- 2 options are available
- Comes with the package of one 5ml bottle and one sponge
- Cost-effective

ANTI-FOG

PRODUCT CODE	DESCRIPTION	PACKAGING
FAF150000	Anti fog solution with alcohol	20 pcs/box
FAF250000	Anti fog solution-alcohol free	20 pcs/box



UNIMAX
UNIMAX MEDICAL SYSTEMS INC.



**UNIMAX
SMOKE
EVACUATION
TUBING**

www.unimaxmeds.com

UNIMAX SMOKE EVACUATION TUBING

SMOKE EVACUATION TUBING

Unimax smoke evacuation tubing offers a simple solution to the removal of smoke caused during laparoscopic procedures. It can not only evacuate smoke efficiently but provide effective protection to the surgical staff.



- 99.999995% efficient of 0.1 micron particles
- A smoke filter and an activated carbon are incorporated
- Lasts throughout a surgical procedure
- Best bacterial and viral control available
- Easy to use

SMOKE EVACUATION TUBING

PRODUCT CODE	DESCRIPTION	PACKAGING
FSE105000	Smoke evacuation tubing	10 pcs/box



UNIMAX
UNIMAX MEDICAL SYSTEMS INC.



UNIMAX MONOPOLAR CABLES

www.unimaxmeds.com

UNIMAX MONOPOLAR CABLES



MONOPOLAR CABLES



8 mm

4 mm

reusable 8 mm

- High quality cable insulation
- Both disposable and reusable cables are available
- Cost-effective

MONOPOLAR CABLES

PRODUCT CODE	DESCRIPTION	PACKAGING
FDC132080	Disposable monopolar cable with 8mm banana connector	40 pcs/box
FDC232080	Disposable monopolar cable with 4mm banana connector	40 pcs/box
TRC132080	Reusable monopolar cable with 8mm banana connector	20 pcs/box



UNIMAX MEDICAL SYSTEMS INC.

8F-2, NO.127 Lane 235, Pao Chiao Rd, Hsin Tien, Taipei, Taiwan
T/+886-2-8919-1698 F/+886-2-8919-1528 www.unimaxmeds.com



UNIMAX
UNIMAX MEDICAL SYSTEMS INC.




UNIMAX ENDOMETRIAL CELL SAMPLER

www.unimaxmeds.com

UNIMAX ENDOMETRIAL CELL SAMPLER



ENDOMETRIAL CELL SAMPLER



Round tip endometrial cell sampler increases patient comfort on insertion

adequate hardness ensures easy insertion

- Allows optimal specimen collection
- Cost-effective

ENDOMETRIAL CELL SAMPLER

PRODUCT CODE	DESCRIPTION	PACKAGING
FEM103000	Endometrial cell sampler 3mm	25 pcs/box
FEM103600	Endometrial cell sampler 3.6mm	25 pcs/box



UNIMAX MEDICAL SYSTEMS INC.

8F-2, NO.127 Lane 235, Pao Chiao Rd, Hsin Tien, Taipei, Taiwan
T/+886-2-8919-1698 F/+886-2-8919-1528 www.unimaxmeds.com



UNIMAX

UNIMAX MEDICAL SYSTEMS INC.



UNIMAX
CLOSED
**SUCTION
CATHETER**

www.unimaxmeds.com

UNIMAX CLOSED SUCTION CATHETER



CLOSED SUCTION CATHETER

- Catheter with centimeter markings allows proper insertion depth in the airway
- Transparent and durable sleeve facilitates catheter manipulation and prevents infiltration of blood and sputum
- ISO color coded button for instant recognition of catheter diameter (ISO 8836: 2009)
- Tethered cap for vacuum connection closing provides full safety
- MDI port provides an effective, quick and convenient way of drug delivery
- Irrigation port with one way valve allows lavage and irrigation which prevents secretion infiltration and ventilation leakage to maintain proper airway pressure
- 360° swivel elbow enables rotation for enhanced patient comfort and reduces torque
- Smooth and soft catheter tip reduces damage to mucous membranes
- Lateral eyes improve airflow and provide diverted fluid path when tip occlusion occurs

Unimax Closed Suction Catheter provides safe and easy access for secretion removal from the airway to minimize infection and maximize diffusion of oxygen through the lung.

CLOSED SUCTION CATHETER

PRODUCT CODE	DESCRIPTION	PACKAGING
FSC105460	Closed suction catheter 5Fr (Grey)	10 pcs/box
FSC106460	Closed suction catheter 6Fr (Light green)	10 pcs/box
FSC107460	Closed suction catheter 7Fr (Ivory)	10 pcs/box
FSC108460	Closed suction catheter 8Fr (Light blue)	10 pcs/box
FSC110460	Closed suction catheter 10Fr (Black)	10 pcs/box
FSC112460	Closed suction catheter 12Fr (White)	10 pcs/box
FSC114460	Closed suction catheter 14Fr (Green)	10 pcs/box
FSC116460	Closed suction catheter 16Fr (Orange)	10 pcs/box




UNIMAX MEDICAL SYSTEMS INC.

8F-2, NO.127 Lane 235, Pao Chiao Rd, Hsin Tien, Taipei, Taiwan
 T/+886-2-8919-1698 F/+886-2-8919-1528 www.unimaxmeds.com



UNIMAX
UNIMAX MEDICAL SYSTEMS INC.



UNIMAX
LIGHT
**HANDLE
COVER**

www.unimaxmeds.com

UNIMAX LIGHT HANDLE COVER



LIGHT HANDLE COVER

Unimax Light handle cover enables easy and quick application which reduces set-up and removal time.



- Smooth material provides universal & secure fit.
- Provides integrated sterile barrier
- Easy to use
- Fits most brands
- Cost effective

LIGHT HANDLE COVER

PRODUCT CODE	DESCRIPTION	PACKAGING
FLH100000	Light handle cover	100 pcs/box



UNIMAX
UNIMAX MEDICAL SYSTEMS INC.



UNIMAX PIPELINE **PROTECTOR**

www.unimaxmeds.com

UNIMAX PIPELINE PROTECTOR



PIPELINE PROTECTOR



Overflow safety filter



Shutoff safety filter



Connector: P type



Connector: D type

- Anti-contamination & bacteria
- Reduces equipment malfunction caused by contamination
- Clear housing provides obvious visual aid for unit replacement
- Comes in 2 filter types: overflow and shutoff safety filter
- 2 types of connector fit most equipments

PIPELINE PROTECTOR

PRODUCT CODE	DESCRIPTION	PACKAGING
1SF100000	Pipeline protector - D type with overflow safety filter	50 pcs/carton
1SF100010	Pipeline protector - D type with shutoff safety filter	50 pcs/carton
1SF200000	Pipeline protector - P type with overflow safety filter	50 pcs/carton
1SF200010	Pipeline protector - P type with shutoff safety filter	50 pcs/carton



UNIMAX


UNIMAX MEDICAL SYSTEMS INC.



UNIMAX SURGICAL MARKING PEN

www.unimaxmeds.com

UNIMAX SURGICAL MARKING PEN



SURGICAL MARKING PEN



Regular point



Fine point

- Scrub resistance
- Non-toxic, non-irritating formula
- Violet ink is visible on most skin tones and prepared skin
- Available with label set and paper ruler
- Comes in 2 options: regular and fine point

SURGICAL MARKING PEN

PRODUCT CODE	DESCRIPTION	PACKAGING
FMP100000	Surgical marking pen – regular point	25 pcs/box
FMP200000	Surgical marking pen – fine point	25 pcs/box